

HSP90 / Heat Shock Protein 90 Mouse anti-Chicken Monoclonal (RPE) (D7A) Antibody - LS-C231930 - LSBio

CatalogID:	LS-C231930
Target:	HSP90 / Heat Shock Protein 90
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	D7A
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	Chicken
Antigen Type:	Purified protein
Immunogen:	Full length protein Hsp90 purified from chicken brain.
Specificity:	Recognizes 90kD proteins corresponding to the molecular mass of Hsp90. Hsp90 alpha specific for human samples. Can isolate complexes of Hsp90, Src kinase and cec37.
Reactivity:	Chicken, Human, Mouse, Rat, Bovine, Pig, Rabbit
Purification:	Protein G purified
Presentation:	PBS, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences