

HSP40 / YDJ1 Mouse anti-Yeast Monoclonal (PerCP) (1G10.H8) Antibody - LS-C231940 - LSBio	
CatalogID:	LS-C231940
Target:	HSP40 / YDJ1
Host	HSP40 / YDJ1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Clone Name:	1G10.H8
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	HSP40 / YDJ1 antibody was raised against Yeast
Antigen Type:	Purified protein
Immunogen:	HSP40 / YDJ1 antibody was raised against full length protein yeast HSP40 (YDJ1).
Specificity:	Detects ~40kD.
Reactivity:	Yeast
Purification:	Protein G purified
Presentation:	0.09% sodium azide, 50% glycerol
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/28/2014 © 2014 LifeSpan BioSciences	