

HSP40 / YDJ1 Mouse anti-Yeast Monoclonal (FITC) (2A7.H6) Antibody - LS-C231954 - LSBio

CatalogID:	LS-C231954
Target:	HSP40 / YDJ1
Host	HSP40 / YDJ1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Clone Name:	2A7.H6
Conjugations:	Fluorescein (FITC)
Immunogen Species:	HSP40 / YDJ1 antibody was raised against Yeast
Antigen Type:	Purified protein
Immunogen:	HSP40 / YDJ1 antibody was raised against full length protein Hsp40 (YDJ1).
Specificity:	Detects ~40kD.
Reactivity:	Yeast
Purification:	Protein G purified
Presentation:	0.09% sodium azide, 50% glycerol
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences