

HSP40 / YDJ1 Mouse anti-Yeast Monoclonal (PerCP) (2A7.H6) Antibody - LS-C231957 - LSBio	
<b>CatalogID:</b>	LS-C231957
<b>Target:</b>	HSP40 / YDJ1
<b>Host</b>	HSP40 / YDJ1 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1,k
<b>Clone Name:</b>	2A7.H6
<b>Conjugations:</b>	Peridinin-chlorophyll-protein (PerCP)
<b>Immunogen Species:</b>	HSP40 / YDJ1 antibody was raised against Yeast
<b>Antigen Type:</b>	Purified protein
<b>Immunogen:</b>	HSP40 / YDJ1 antibody was raised against full length protein Hsp40 (YDJ1).
<b>Specificity:</b>	Detects ~40kD.
<b>Reactivity:</b>	Yeast
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	0.09% sodium azide, 50% glycerol
<b>Recommended Storage:</b>	Store at -20°C.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	