

Rhodopsin / RHO Mouse anti-Bovine Monoclonal (PerCP) (4D2) Antibody - LS-C231973 - LSBio	
CatalogID:	LS-C231973
Target:	rhodopsin
Synonyms:	RHO Antibody, CSNBAD1 Antibody, OPN2 Antibody, Rhodopsin Antibody, RP4 Antibody, Opsin 2, rod pigment Antibody, Opsin-2 Antibody
Family / Subfamily:	GPCR / Opsin
Host	RHO antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	4D2
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	Rhodopsin / RHO antibody was raised against Bovine
Antigen Type:	Purified protein
Immunogen:	Rhodopsin / RHO antibody was raised against bovine rhodopsin.
Specificity:	N-terminus of Rhodopsin.
Reactivity:	Bovine, Mouse, Rat, Pig, Xenopus, Fish, Bird
Non-Reactivity:	Drosophila, Invertebrate
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014 © 2014 LifeSpan BioSciences	