

HSPA5 / GRP78 / BIP Mouse anti-Human Monoclonal (PerCP) (1H11-1H7) Antibody - LS-C232042 - LSBio

CatalogID:	LS-C232042
Target:	heat shock 70kDa protein 5 (glucose-regulated protein, 78kDa) (HSPA5)
Synonyms:	HSPA5 Antibody, BIP Antibody, Heat shock 70 kDa protein 5 Antibody, MIF2 Antibody, GRP-78 Antibody, GRP78 Antibody
Host	HSPA5 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	1H11-1H7
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	HSPA5 / GRP78 / BIP antibody was raised against Human
Antigen Type:	Fusion protein
Immunogen:	HSPA5 / GRP78 / BIP antibody was raised against his-tagged human Grp78
Specificity:	Detects a protein with a mass of ~78kD corresponding to Grp78.
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences