

RAB5 Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C232073 - LSBio

CatalogID:	LS-C232073
Target:	Ras-related protein Rab-5 (RAB5)
Host	Ras-related protein Rab-5 (RAB5) antibody was produced in Rabbit
Clonality:	Polyclonal
Conjugations:	Fluorescein (FITC)
Immunogen Species:	RAB5 antibody was raised against Human
Antigen Type:	Synthetic peptide - KLH conjugated
Immunogen:	RAB5 antibody was raised against human Rab5 synthetic peptide conjugated to KLH; identical to dog Rab5 sequence over the residues.
Specificity:	The antibody recognizes 26kD Rab5 of the above species origins.
Reactivity:	Human, Monkey, Mouse, Rat, Bovine
Purification:	Affinity purified
Presentation:	PBS, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Immunofluorescence, Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences