

RAB5 Rabbit anti-Human Polyclonal (PerCP) Antibody - LS-C232076 - LSBio	
<b>CatalogID:</b>	LS-C232076
<b>Target:</b>	Ras-related protein Rab-5 (RAB5)
<b>Host</b>	Ras-related protein Rab-5 (RAB5) antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Conjugations:</b>	Peridinin-chlorophyll-protein (PerCP)
<b>Immunogen Species:</b>	RAB5 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide - KLH conjugated
<b>Immunogen:</b>	RAB5 antibody was raised against human Rab5 synthetic peptide conjugated to KLH; identical to dog Rab5 sequence over the residues.
<b>Specificity:</b>	The antibody recognizes 26kD Rab5 of the above species origins.
<b>Reactivity:</b>	Human, Monkey, Mouse, Rat, Bovine
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, 0.09% sodium azide, 50% glycerol.
<b>Recommended Storage:</b>	Store at -20°C.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Immunofluorescence, Western blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	