

CANX / Calnexin Rabbit anti-Human Polyclonal (C-Terminus) (PerCP) Antibody - LS-C232111 - LSBio	
CatalogID:	LS-C232111
Target:	calnexin (CANX)
Synonyms:	CANX Antibody, CNX Antibody, Calnexin Antibody, IP90 Antibody, p90 Antibody
Host	CANX antibody was produced in Rabbit
Clonality:	Polyclonal
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	CANX / Calnexin antibody was raised against Human
Antigen Type:	Synthetic peptide - KLH conjugated
Immunogen:	CANX / Calnexin antibody was raised against dog Calnexin C-terminal synthetic peptide conjugated to KLH. Identical to human, mouse and rat calnexin sequences over these residues.
Specificity:	Detects ~90kD.
Epitope:	C-Terminus
Reactivity:	Human, Monkey, Mouse, Rat, Bovine, Dog, Guinea pig, Hamster, Pig, Rabbit, Sheep, Chicken, Xenopus, Drosophila, Quail
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	200 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/28/2014 © 2014 LifeSpan BioSciences	