

HSPE1 / HSP10 / Chaperonin 10 Rabbit anti-Human Polyclonal (aa91-101) (PerCP) Antibody - LS-C232127 - LSBio	
CatalogID:	LS-C232127
Target:	heat shock 10kDa protein 1 (HSPE1)
Synonyms:	HSPE1 Antibody, 10 kd chaperonin Antibody, 10 kDa chaperonin Antibody, Chaperonin 10 Antibody, CPN10 Antibody, Early-pregnancy factor Antibody, EPF Antibody, GROES Antibody, HSP10 Antibody
Host	HSPE1 antibody was produced in Rabbit
Clonality:	Polyclonal
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	HSPE1 / HSP10 / Chaperonin 10 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	HSPE1 / HSP10 / Chaperonin 10 antibody was raised against human Cpn10 peptide AA 91-101.
Specificity:	Detects ~10kD.
Epitope:	aa91-101
Reactivity:	Human, Mouse, Rat, Bovine, Dog, Guinea pig, Pig, Rabbit, Sheep, Xenopus
Purification:	Protein A purified
Presentation:	PBS, pH 7.3, 0.02% sodium azide.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Immunofluorescence, Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	