

HSP70 / Heat Shock Protein 70 Chicken anti-Human Polyclonal (FITC) Antibody - LS-C232141 - LSBio

CatalogID:	LS-C232141
Target:	HSP70 / Heat Shock Protein 70
Host	Chicken
Clonality:	Polyclonal
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Human
Antigen Type:	Purified protein
Immunogen:	Full length Human Hsp70.
Specificity:	Hsp70.
Reactivity:	Human, Mouse, Rat
Purification:	PEG purified
Presentation:	PBS, pH 7.4, 50% glycerol, 0.09% sodium azide
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences