

| HSP70 / Heat Shock Protein 70 Chicken anti-Human Polyclonal (PerCP) Antibody - LS-C232144 - LSBio | |
|--|--|
| CatalogID: | LS-C232144 |
| Target: | HSP70 / Heat Shock Protein 70 |
| Host | Chicken |
| Clonality: | Polyclonal |
| Conjugations: | Peridinin-chlorophyll-protein (PerCP) |
| Immunogen Species: | Human |
| Antigen Type: | Purified protein |
| Immunogen: | Full length Human Hsp70. |
| Specificity: | Hsp70. |
| Reactivity: | Human, Mouse, Rat |
| Purification: | PEG purified |
| Presentation: | PBS, pH 7.4, 50% glycerol, 0.09% sodium azide |
| Recommended Storage: | Store at -20°C. |
| Usage Summary: | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested. |
| Uses: | Western blot, ELISA (Optimal dilution to be determined by the researcher) |
| Size: | 100 µg |
| Requested From: | Japan |
| Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/28/2014 © 2014 LifeSpan BioSciences | |