

**HSP70 / Heat Shock Protein 70 Chicken anti-Human Polyclonal (RPE) Antibody - LS-C232153 - LSBio**

<b>CatalogID:</b>	LS-C232153
<b>Target:</b>	HSP70 / Heat Shock Protein 70
<b>Host</b>	Chicken
<b>Clonality:</b>	Polyclonal
<b>Conjugations:</b>	R. Phycoerythrin (RPE)
<b>Immunogen Species:</b>	Human
<b>Antigen Type:</b>	Purified protein
<b>Immunogen:</b>	Full length Human Hsp70.
<b>Specificity:</b>	Hsp70.
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Purification:</b>	PEG purified
<b>Presentation:</b>	PBS, pH 7.4, 50% glycerol, 0.09% sodium azide
<b>Recommended Storage:</b>	Store at -20°C.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences