

**AP1S1 Rabbit anti-Human Polyclonal (aa15-45) (APC) Antibody - LS-C232716 - LSBio**

|                             |  |
|-----------------------------|--|
| <b>CatalogID:</b>           | LS-C232716   |
| <b>Target:</b>              | adaptor-related protein complex 1, sigma 1 subunit (AP1S1)   |
| <b>Synonyms:</b>            | AP1S1 Antibody, AP19 Antibody, AP-1 complex subunit sigma-1A Antibody, CLAPS1 Antibody, HA1 19 kDa subunit Antibody, SIGMA1A Antibody, Sigma1A-adaptin Antibody, WUGSC:H_DJ0747G18.2 Antibody, Sigma-adaptin 1A Antibody |
| <b>Host</b>                 | AP1S1 antibody was produced in Rabbit  |
| <b>Clonality:</b>           | Polyclonal   |
| <b>Isotype:</b>             | IgG  |
| <b>Conjugations:</b>        | Allophycocyanin (APC)  |
| <b>Immunogen Species:</b>   | AP1S1 antibody was raised against Human  |
| <b>Antigen Type:</b>        | Synthetic peptide  |
| <b>Immunogen:</b>           | AP1S1 antibody was raised against kLH-conjugated synthetic peptide (aa15-45) from the N-terminal region of human AP1S1.  |
| <b>Specificity:</b>         | Human AP1S1.   |
| <b>Epitope:</b>             | aa15-45  |
| <b>Reactivity:</b>          | Human  |
| <b>Purification:</b>        | Affinity purified  |
| <b>Presentation:</b>        | PBS, pH 7.2, 0.09% sodium azide.   |
| <b>Recommended Storage:</b> | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.  |
| <b>Usage Summary:</b>       | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.   |
| <b>Uses:</b>                | IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)   |
| <b>Size:</b>                | 200 µl   |
| <b>Requested From:</b>      | Japan  |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences