

AP1S1 Rabbit anti-Human Polyclonal (aa15-45) (AP) Antibody - LS-C232717 - LSBio

CatalogID:	LS-C232717
Target:	adaptor-related protein complex 1, sigma 1 subunit (AP1S1)
Synonyms:	AP1S1 Antibody, AP19 Antibody, AP-1 complex subunit sigma-1A Antibody, CLAPS1 Antibody, HA1 19 kDa subunit Antibody, SIGMA1A Antibody, Sigma1A-adaptin Antibody, WUGSC:H_DJ0747G18.2 Antibody, Sigma-adaptin 1A Antibody
Host	AP1S1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Alkaline Phosphatase (AP)
Immunogen Species:	AP1S1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	AP1S1 antibody was raised against kLH-conjugated synthetic peptide (aa15-45) from the N-terminal region of human AP1S1.
Specificity:	Human AP1S1.
Epitope:	aa15-45
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences