

AP3M1 Rabbit anti-Human Polyclonal (aa389-418) (AP) Antibody - LS-C232722 - LSBio

CatalogID:	LS-C232722
Target:	Adaptor-related Protein Complex 3, Mu 1 Subunit (AP3M1)
Synonyms:	AP3M1 Antibody, AP-3 complex subunit mu-1 Antibody, Mu3A-adaptin Antibody, Mu-adaptin 3A Antibody
Host	AP3M1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Alkaline Phosphatase (AP)
Immunogen Species:	AP3M1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	AP3M1 antibody was raised against kLH-conjugated synthetic peptide (aa389-418) from the C-terminus of human AP3M1.
Specificity:	Human AP3M1
Epitope:	aa389-418
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences