

**AP4M1 Rabbit anti-Human Polyclonal (aa222-251) (APC) Antibody - LS-C232726 - LSBio**

<b>CatalogID:</b>	LS-C232726
<b>Target:</b>	adaptor-related protein complex 4, mu 1 subunit (AP4M1)
<b>Synonyms:</b>	AP4M1 Antibody, AP-4 complex subunit mu-1 Antibody, MU-4 Antibody, MU-ARP2 Antibody, Mu-adaptin-related protein 2 Antibody, Mu4-adaptin Antibody, MUARP2 Antibody, SPG50 Antibody, Mu subunit of AP-4 Antibody, Mu-adaptin-related protein-2 Antibody, Mu4 Antibody, CPSQ3 Antibody
<b>Host</b>	AP4M1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Allophycocyanin (APC)
<b>Immunogen Species:</b>	AP4M1 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	AP4M1 antibody was raised against kLH-conjugated synthetic peptide (aa222-251) from internal region of human AP4M1.
<b>Specificity:</b>	Human AP4M1.
<b>Epitope:</b>	aa222-251
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences