

AP4M1 Rabbit anti-Human Polyclonal (aa222-251) (AP) Antibody - LS-C232729 - LSBio

CatalogID:	LS-C232729
Target:	adaptor-related protein complex 4, mu 1 subunit (AP4M1)
Synonyms:	AP4M1 Antibody, AP-4 complex subunit mu-1 Antibody, MU-4 Antibody, MU-ARP2 Antibody, Mu-adaptin-related protein 2 Antibody, Mu4-adaptin Antibody, MUARP2 Antibody, SPG50 Antibody, Mu subunit of AP-4 Antibody, Mu-adaptin-related protein-2 Antibody, Mu4 Antibody, CPSQ3 Antibody
Host	AP4M1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Alkaline Phosphatase (AP)
Immunogen Species:	AP4M1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	AP4M1 antibody was raised against kLH-conjugated synthetic peptide (aa222-251) from internal region of human AP4M1.
Specificity:	Human AP4M1.
Epitope:	aa222-251
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences