

ASAP3 Rabbit anti-Human Polyclonal (aa26-54) (PE) Antibody - LS-C232889 - LSBio

CatalogID:	LS-C232889
Target:	ArfGAP with SH3 domain, ankyrin repeat and PH domain 3 (ASAP3)
Synonyms:	ASAP3 Antibody, ACAP4 Antibody, ARF6 GTPase-activating protein Antibody, CENTB6 Antibody, DDEFL1 Antibody, Centaurin, beta 6 Antibody, UPLC1 Antibody
Host	ASAP3 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ASAP3 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ASAP3 antibody was raised against kLH-conjugated synthetic peptide (aa26-54) from the N-terminal region of human ASAP3.
Specificity:	Human ASAP3
Epitope:	aa26-54
Reactivity:	Human, Hamster
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences