

ATP11C Rabbit anti-Human Polyclonal (aa596-625) (FITC) Antibody - LS-C232980 - LSBio	
CatalogID:	LS-C232980
Target:	ATPase, class VI, type 11C (ATP11C)
Synonyms:	ATP11C Antibody, ATPase IQ Antibody, ATPase, class VI, type 11C Antibody, ATPIG Antibody, ATPase class VI type 11C Antibody, ATPIQ Antibody
Family / Subfamily:	Transporter / ATPase - P type, type IV
Host	ATP11C antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ATP11C antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ATP11C antibody was raised against kLH-conjugated synthetic peptide (aa596-625) from internal region of human ATP11C.
Specificity:	Human ATP11C
Epitope:	aa596-625
Reactivity:	Human, Mouse, Hamster
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	