

BMPER Rabbit anti-Human Polyclonal (aa155-186) (FITC) Antibody - LS-C233180 - LSBio

CatalogID:	LS-C233180
Target:	BMP binding endothelial regulator (BMPER)
Synonyms:	BMPER Antibody, CRIM3 Antibody, Crossveinless 2 Antibody, CV-2 Antibody, Crossveinless-2 Antibody, HCV2 Antibody, Protein crossveinless-2 Antibody, CV2 Antibody, KIAA1965 Antibody
Host	BMPER antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	BMPER antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	BMPER antibody was raised against kLH-conjugated synthetic peptide (aa155-186) from the N-terminal region of human BMPER.
Specificity:	Human BMPER
Epitope:	aa155-186
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences