

CCNB3 / Cyclin B3 Rabbit anti-Rat Polyclonal (pSer284) (PE) Antibody - LS-C233580 - LSBio	
CatalogID:	LS-C233580
Target:	cyclin B3 (CCNB3)
Synonyms:	CCNB3 Antibody, Cyclin B3 Antibody, G2/mitotic-specific cyclin-B3 Antibody, CYCB3 Antibody
Host	CCNB3 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CCNB3 / Cyclin B3 antibody was raised against Rat
Antigen Type:	Synthetic phospho-peptide
Immunogen:	CCNB3 / Cyclin B3 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser284 of rat CCNB3.
Specificity:	Rat CCNB3 / Cyclin B3.
Epitope:	pSer284
Reactivity:	Rat, Rabbit
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	