

CCNYL1 Rabbit anti-Human Polyclonal (N-Terminus) (PE) Antibody - LS-C233627 - LSBio

CatalogID:	LS-C233627
Target:	cyclin Y-like 1 (CCNYL1)
Synonyms:	CCNYL1 Antibody, Cyclin Y-like 1 Antibody, Cyclin-Y-like protein 1 Antibody
Host	CCNYL1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CCNYL1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CCNYL1 antibody was raised against kLH-conjugated synthetic peptide between 48-78 amino acids from the N-terminal region of human CCYL1.
Specificity:	Human CCNYL1
Epitope:	N-Terminus
Reactivity:	Human, Hamster
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences