

CREB1 / CREB Rabbit anti-Human Polyclonal (aa112-140) (FITC) Antibody - LS-C234220 - LSBio	
CatalogID:	LS-C234220
Target:	cAMP responsive element binding protein 1 (CREB1)
Synonyms:	CREB1 Antibody, CREB-1 Antibody, Transactivator protein Antibody, CREB Antibody
Family / Subfamily:	Transcription factor / not assigned-Transcription factor
Host	CREB1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CREB1 / CREB antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CREB1 / CREB antibody was raised against kLH-conjugated synthetic peptide (aa112-140) from internal region of human CREB1.
Specificity:	Human CREB1 / CREB
Epitope:	aa112-140
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ICC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	