

CD95 / FAS Mouse anti-Human Monoclonal (FITC) Antibody - LS-C23440 - LSBio

CatalogID:	LS-C23440
Target:	Fas cell surface death receptor
Synonyms:	FAS Antibody, ALPS1A Antibody, APO-1 Antibody, APT1 Antibody, Apo-1 antigen Antibody, APO-1 cell surface antigen Antibody, Apo11 antigen Antibody, Apoptosis antigen 1 Antibody, CD95 antigen Antibody, Delta Fas/APO-1/CD95 Antibody, Fas AMA Antibody, FAS1 Antibody, FasI receptors Antibody, FASLG receptor Antibody, FASTM Antibody, FAS 827dupA Antibody, FasI receptor Antibody, TNFRSF6 Antibody, CD95 Antibody
Family / Subfamily:	TNF Receptor / not assigned-TNF Receptor
Host	FAS antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD95 / FAS antibody was raised against Human
Specificity:	Recognizes human Fas (CD95, APO-1).
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, sodium azide.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Immunohistochemistry (frozen sections). Immunoprecipitation. In vitro induction of apoptosis. Flow Cytometry: 10 ul/1E6 cells.
Uses:	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences