

CD95 / FAS Mouse anti-Human Monoclonal (PE) Antibody - LS-C23448 - LSBio

<b>CatalogID:</b>	LS-C23448
<b>Target:</b>	Fas cell surface death receptor
<b>Synonyms:</b>	FAS Antibody, ALPS1A Antibody, APO-1 Antibody, APT1 Antibody, Apo-1 antigen Antibody, APO-1 cell surface antigen Antibody, Apo11 antigen Antibody, Apoptosis antigen 1 Antibody, CD95 antigen Antibody, Delta Fas/APO-1/CD95 Antibody, Fas AMA Antibody, FAS1 Antibody, FasI receptors Antibody, FASLG receptor Antibody, FASTM Antibody, FAS 827dupA Antibody, FasI receptor Antibody, TNFRSF6 Antibody, CD95 Antibody
<b>Family / Subfamily:</b>	TNF Receptor / not assigned-TNF Receptor
<b>Host</b>	FAS antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	CD95 / FAS antibody was raised against Human
<b>Specificity:</b>	Recognizes human Fas (CD95, APO-1).
<b>Reactivity:</b>	Human
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, 0.1% sodium azide, stabilizing agent.
<b>Recommended Storage:</b>	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
<b>Usage Summary:</b>	Flow cytometry. Immunohistochemistry (frozen sections). Immunoprecipitation. In vitro induction of apoptosis. Flow Cytometry: 10 ul/1E6 cells.
<b>Uses:</b>	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 tst
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences