

ERAS Rabbit anti-Human Polyclonal (aa20-50) (PE) Antibody - LS-C235011 - LSBio

CatalogID:	LS-C235011
Target:	ES cell expressed Ras (ERAS)
Synonyms:	ERAS Antibody, ES cell expressed Ras Antibody, E-Ras Antibody, GTPase ERas Antibody, HRASP Antibody, HRAS2 Antibody, Small GTPase protein E-Ras Antibody
Host	ERAS antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ERAS antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ERAS antibody was raised against kLH-conjugated synthetic peptide (aa20-50) from the N-terminal region of human ERAS.
Specificity:	Human ERAS
Epitope:	aa20-50
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences