

FBXO34 Rabbit anti-Human Polyclonal (C-Terminus) (PE) Antibody - LS-C235232 - LSBio

CatalogID:	LS-C235232
Target:	F-box protein 34 (FBXO34)
Synonyms:	FBXO34 Antibody, CGI-301 Antibody, F-box only protein 34 Antibody, Fbx34 Antibody, F-box protein 34 Antibody, Protein CGI-301 Antibody
Family / Subfamily:	F-Box
Host	FBXO34 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	FBXO34 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	FBXO34 antibody was raised against kLH-conjugated synthetic peptide between 607-637 amino acids from the C-terminus of human FBXO34.
Specificity:	Human FBXO34
Epitope:	C-Terminus
Reactivity:	Human, Hamster
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences