

**GNAS Rabbit anti-Human Polyclonal (aa286-315) (PE) Antibody - LS-C235787 - LSBio**

<b>CatalogID:</b>	LS-C235787
<b>Target:</b>	GNAS complex locus
<b>Synonyms:</b>	GNAS Antibody, AHO Antibody, C20orf45 Antibody, Extra large alphas protein Antibody, GSP Antibody, GNAS complex locus Antibody, GNASXL Antibody, NESP Antibody, PHP1B Antibody, PHP1C Antibody, PHP1A Antibody, POH Antibody, SCG6 Antibody, XLalphas Antibody, GNAS1 Antibody, GPSA Antibody, GSA Antibody, NESP55 Antibody, Protein ALEX Antibody, Secretogranin VI Antibody
<b>Host</b>	GNAS antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	GNAS antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	GNAS antibody was raised against kLH-conjugated synthetic peptide between 286 -315aa from the C-terminus of human GNAS.
<b>Specificity:</b>	Human GNAS
<b>Epitope:</b>	aa286-315
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, ICC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences