

GPAM Rabbit anti-Human Polyclonal (aa433-463) (FITC) Antibody - LS-C235870 - LSBio

CatalogID:	LS-C235870
Target:	glycerol-3-phosphate acyltransferase, mitochondrial (GPAM)
Synonyms:	GPAM Antibody, GPAT Antibody, GPAT1 Antibody, GPAT-1 Antibody, KIAA1560 Antibody, RP11-426E5.2 Antibody
Host	GPAM antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	GPAM antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	GPAM antibody was raised against kLH-conjugated synthetic peptide (aa433-463) from internal region of human GPAM.
Specificity:	Human GPAM
Epitope:	aa433-463
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences