

GRINA Rabbit anti-Human Polyclonal (aa130-158) (PE) Antibody - LS-C235998 - LSBio

CatalogID:	LS-C235998
Target:	glutamate receptor, ionotropic, N-methyl D-aspartate-associated protein 1 (glutamate binding) (GRINA)
Synonyms:	GRINA Antibody, Protein lifeguard 1 Antibody, LFG1 Antibody, TMBIM3 Antibody, HNRGW Antibody, NMDARA1 Antibody
Host	GRINA antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	GRINA antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	GRINA antibody was raised against kLH-conjugated synthetic peptide (aa130-158) from internal region of human GRINA.
Specificity:	Human GRINA
Epitope:	aa130-158
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences