

GTPBP2 Rabbit anti-Human Polyclonal (aa543-572) (FITC) Antibody - LS-C236044 - LSBio

CatalogID:	LS-C236044
Target:	GTP binding protein 2 (GTPBP2)
Synonyms:	GTPBP2 Antibody, GTP binding protein 2 Antibody, GTP-binding protein 2 Antibody
Host	GTPBP2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	GTPBP2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	GTPBP2 antibody was raised against kLH-conjugated synthetic peptide (aa543-572) from the C-terminus of human GTPBP2.
Specificity:	Human GTPBP2
Epitope:	aa543-572
Reactivity:	Human, Mouse, Hamster
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences