

HP / Haptoglobin Rabbit anti-Human Polyclonal (aa299-326) (PE) Antibody - LS-C236390 - LSBio	
CatalogID:	LS-C236390
Target:	haptoglobin (HP)
Synonyms:	HP Antibody, Binding peptide Antibody, Haptoglobin alpha(1S)-beta Antibody, Haptoglobin, beta polypeptide Antibody, Haptoglobin Antibody, Haptoglobin alpha(2FS)-beta Antibody, Haptoglobin, alpha polypeptide Antibody, HP2ALPHA2 Antibody, BP Antibody, Zonulin Antibody, HPA1S Antibody
Host	HP antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	HP / Haptoglobin antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	HP / Haptoglobin antibody was raised against kLH-conjugated synthetic peptide (aa299-326) from internal region of human HP.
Specificity:	Human HP / Haptoglobin
Epitope:	aa299-326
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	