

HPR Rabbit anti-Human Polyclonal (aa219-247) (PE) Antibody - LS-C236392 - LSBio

CatalogID:	LS-C236392
Target:	Haptoglobin-related Protein (HPR)
Synonyms:	HPR Antibody, A-259H10.2 Antibody, Haptoglobin-related protein Antibody, Haptoglobin-related locus Antibody
Host	HPR antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	HPR antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	HPR antibody was raised against kLH-conjugated synthetic peptide (aa219-247) from internal region of human HPR.
Specificity:	Human HPR
Epitope:	aa219-247
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences