

MYCT1 Rabbit anti-Human Polyclonal (aa142-171) (Biotin) Antibody - LS-C237860 - LSBio	
CatalogID:	LS-C237860
Target:	Myc Target 1 (MYCT1)
Synonyms:	MYCT1 Antibody, MTMC1 Antibody, Myc target protein 1 Antibody, Myc target 1 Antibody, Myc target in myeloid cells 1 Antibody, MTLC Antibody
Host	MYCT1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Biotin
Immunogen Species:	MYCT1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MYCT1 antibody was raised against kLH-conjugated synthetic peptide (aa142-171) from the C-terminus of human MYCT1.
Specificity:	Human MYCT1
Epitope:	aa142-171
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	