

NEU2 Rabbit anti-Human Polyclonal (aa30-58) (FITC) Antibody - LS-C237971 - LSBio

CatalogID:	LS-C237971
Target:	sialidase 2 (cytosolic sialidase) (NEU2)
Synonyms:	NEU2 Antibody, N-acetyl-alpha-neuraminidase 2 Antibody, Sialidase-2 Antibody, SIAL2 Antibody, Cytosolic sialidase Antibody, Neuraminidase 2 Antibody
Host	NEU2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	NEU2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	NEU2 antibody was raised against kLH-conjugated synthetic peptide (aa30-58) from the N-terminal region of human NEU2.
Specificity:	Human NEU2
Epitope:	aa30-58
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences