

PARVA Rabbit anti-Mouse Polyclonal (aa65-93) (FITC) Antibody - LS-C238700 - LSBio

CatalogID:	LS-C238700
Target:	parvin, alpha (PARVA)
Synonyms:	PARVA Antibody, Actopaxin Antibody, Alpha-parvin Antibody, CH-ILKBP Antibody, Parvin, alpha Antibody, MXRA2 Antibody
Host	PARVA antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PARVA antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	PARVA antibody was raised against kLH-conjugated synthetic peptide (aa65-93) from the N-terminal region of mouse Parva.
Specificity:	Mouse PARVA.
Epitope:	aa65-93
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences