

PCDHB3 Rabbit anti-Human Polyclonal (aa19-47) (PE) Antibody - LS-C238817 - LSBio

CatalogID:	LS-C238817
Target:	Protocadherin Beta 3 (PCDHB3)
Synonyms:	PCDHB3 Antibody, Protocadherin beta-3 Antibody, Protocadherin beta 3 Antibody, PCDH-beta-3 Antibody, PCDH-BETA3 Antibody
Family / Subfamily:	Cadherin / not assigned-Cadherin
Host	PCDHB3 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	PCDHB3 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PCDHB3 antibody was raised against kLH-conjugated synthetic peptide (aa19-47) from the N-terminal region of human PCDHB3.
Specificity:	Human PCDHB3.
Epitope:	aa19-47
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences