

PDC Rabbit anti-Human Polyclonal (aa58-89) (FITC) Antibody - LS-C238870 - LSBio

CatalogID:	LS-C238870
Target:	Phosducin (PDC)
Synonyms:	PDC Antibody, G beta gamma binding protein Antibody, MEKA Antibody, PHD Antibody, Phosducin-like orphan protein Antibody, PhLP Antibody, Phosducin Antibody, PhLOP Antibody, Protein MEKA Antibody
Host	PDC antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PDC antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PDC antibody was raised against kLH-conjugated synthetic peptide (aa58-89) from internal region of human PDC.
Specificity:	Human PDC
Epitope:	aa58-89
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences