

PERP Rabbit anti-Human Polyclonal (aa161-191) (FITC) Antibody - LS-C238924 - LSBio

CatalogID:	LS-C238924
Target:	PERP, TP53 apoptosis effector
Synonyms:	PERP Antibody, DJ496H19.1 Antibody, KCP-1 Antibody, PERP, TP53 apoptosis effector Antibody, KCP1 Antibody, p53-induced protein PIGPC1 Antibody, THW Antibody, KRTCAP1 Antibody, PIGPC1 Antibody, RP3-496H19.1 Antibody, Transmembrane protein THW Antibody
Host	PERP antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PERP antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PERP antibody was raised against kLH-conjugated synthetic peptide (aa161-191) from the C-terminus of human PERP.
Specificity:	Human PERP.
Epitope:	aa161-191
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences