

PREPL Rabbit anti-Human Polyclonal (aa631-660) (PE) Antibody - LS-C239427 - LSBio

CatalogID:	LS-C239427
Target:	Prolyl Endopeptidase-like (PREPL)
Synonyms:	PREPL Antibody, KIAA0436 Antibody, Prolyl endopeptidase-like Antibody, Prolylendopeptidase-like Antibody, Putative prolyl oligopeptidase Antibody
Family / Subfamily:	Protease / Serine S9
Host	PREPL antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	PREPL antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PREPL antibody was raised against kLH-conjugated synthetic peptide (aa631-660) from the C-terminus of human PREPL.
Specificity:	Human PREPL
Epitope:	aa631-660
Reactivity:	Human, Mouse, Hamster
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences