

PTPLA Rabbit anti-Human Polyclonal (aa36-66) (FITC) Antibody - LS-C239595 - LSBio

CatalogID:	LS-C239595
Target:	protein tyrosine phosphatase-like (proline instead of catalytic arginine), member A (PTPLA)
Synonyms:	PTPLA Antibody, CAP Antibody, Cementum attachment protein Antibody, HACD1 Antibody
Family / Subfamily:	Protein Phosphatase / PTP-like
Host	PTPLA antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PTPLA antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PTPLA antibody was raised against KLH-conjugated synthetic peptide (aa36-66) from the N-terminal region of human PTPLA.
Specificity:	Human PTPLA
Epitope:	aa36-66
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences