

RBPJL Rabbit anti-Human Polyclonal (aa15-44) (PE) Antibody - LS-C239856 - LSBio

CatalogID:	LS-C239856
Target:	recombination signal binding protein for immunoglobulin kappa J region-like (RBPJL)
Synonyms:	RBPJL Antibody, RBP-L Antibody, RBPSUHL Antibody, SUHL Antibody, SUH Antibody, RBPL Antibody, Transcription factor RBP-L Antibody
Host	RBPJL antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	RBPJL antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	RBPJL antibody was raised against kLH-conjugated synthetic peptide (aa15-44) from the N-terminal region of human RBPJL.
Specificity:	Human RBPJL
Epitope:	aa15-44
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences