

RET Rabbit anti-Human Polyclonal (aa159-189) (PE) Antibody - LS-C239931 - LSBio

CatalogID:	LS-C239931
Target:	ret proto-oncogene
Synonyms:	RET Antibody, Cadherin family member 12 Antibody, CDHF12 Antibody, C-ret Antibody, Hirschsprung disease 1 Antibody, HSCR1 Antibody, Hydroxyaryl-protein kinase Antibody, Oncogene ret Antibody, Proto-oncogene c-Ret Antibody, RET-ELE1 Antibody, RET51 Antibody, Ret proto-oncogene Antibody, RET transforming sequence Antibody, RET43 Antibody, Men2-ret Antibody, MEN2B Antibody, MTC1 Antibody, Receptor tyrosine kinase Antibody, CDHR16 Antibody, MEN2A Antibody, PTC Antibody, RET9 Antibody
Family / Subfamily:	Protein Kinase / Cadherin
Host	RET antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	RET antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	RET antibody was raised against kLH-conjugated synthetic peptide (aa159-189) from the N-terminal region of human RET.
Specificity:	Human RET
Epitope:	aa159-189
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences