

RHOC Rabbit anti-Human Polyclonal (aa101-130) (PE) Antibody - LS-C239977 - LSBio

CatalogID:	LS-C239977
Target:	ras homolog family member C (RHOC)
Synonyms:	RHOC Antibody, ARHC Antibody, Rho cDNA clone 9 Antibody, RhoC GTPase Antibody, RHOH9 Antibody, Small GTP binding protein RhoC Antibody, ARH9 Antibody, H9 Antibody, Oncogene RHO H9 Antibody, Ras homolog family member C Antibody, RAS-related homolog 9 Antibody
Host	RHOC antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	RHOC antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	RHOC antibody was raised against kLH-conjugated synthetic peptide (aa101-130) from the C-terminus of human RHOC.
Specificity:	Human RHOC
Epitope:	aa101-130
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences