

RILPL2 Rabbit anti-Human Polyclonal (aa120-149) (FITC) Antibody - LS-C240020 - LSBio	
CatalogID:	LS-C240020
Target:	Rab Interacting Lysosomal Protein-Like 2 (RILPL2)
Synonyms:	RILPL2 Antibody, p40phox-binding protein Antibody, RLP2 Antibody, RILP-like protein 2 Antibody
Host	RILPL2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	RILPL2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	RILPL2 antibody was raised against kLH-conjugated synthetic peptide (aa120-149) from internal region of human RILPL2.
Specificity:	Human RILPL2.
Epitope:	aa120-149
Reactivity:	Human, Hamster
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	