

TWISTNB Rabbit anti-Human Polyclonal (aa306-336) (PE) Antibody - LS-C240191 - LSBio	
CatalogID:	LS-C240191
Target:	TWIST neighbor (TWISTNB)
Synonyms:	TWISTNB Antibody, TWIST neighbor Antibody, Twist neighbor protein Antibody
Host	TWISTNB antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	TWISTNB antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	TWISTNB antibody was raised against kLH-conjugated synthetic peptide (aa306-336) from the C-terminus of human RPA43.
Specificity:	Human TWISTNB
Epitope:	aa306-336
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/28/2014 © 2014 LifeSpan BioSciences	