

SCEL Rabbit anti-Human Polyclonal (aa247-277) (FITC) Antibody - LS-C240438 - LSBio

CatalogID:	LS-C240438
Target:	sciellin (SCEL)
Synonyms:	SCEL Antibody, Sciellin Antibody
Host	SCEL antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	SCEL antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	SCEL antibody was raised against kLH-conjugated synthetic peptide (aa247-277) from internal region of human SCEL.
Specificity:	Human SCEL
Epitope:	aa247-277
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences