

**SOCS4 Rabbit anti-Human Polyclonal (aa234-262) (FITC) Antibody - LS-C241108 - LSBio**

|                             |   |
|-----------------------------|---|
| <b>CatalogID:</b>           | LS-C241108  |
| <b>Target:</b>              | suppressor of cytokine signaling 4 (SOCS4)  |
| <b>Synonyms:</b>            | SOCS4 Antibody, CIS4 Antibody, SOCS-4 Antibody, SSI4 Antibody, SOCS-7 Antibody  |
| <b>Host</b>                 | SOCS4 antibody was produced in Rabbit   |
| <b>Clonality:</b>           | Polyclonal  |
| <b>Isotype:</b>             | IgG   |
| <b>Conjugations:</b>        | Fluorescein (FITC)  |
| <b>Immunogen Species:</b>   | SOCS4 antibody was raised against Human   |
| <b>Antigen Type:</b>        | Synthetic peptide   |
| <b>Immunogen:</b>           | SOCS4 antibody was raised against kLH-conjugated synthetic peptide (aa234-262) from internal region of human SOCS4.                         |
| <b>Specificity:</b>         | Human SOCS4.  |
| <b>Epitope:</b>             | aa234-262   |
| <b>Reactivity:</b>          | Human   |
| <b>Purification:</b>        | Affinity purified   |
| <b>Presentation:</b>        | PBS, pH 7.2, 0.09% sodium azide.  |
| <b>Recommended Storage:</b> | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months. |
| <b>Usage Summary:</b>       | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.                  |
| <b>Uses:</b>                | IHC, ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)   |
| <b>Size:</b>                | 200 µl  |
| <b>Requested From:</b>      | Japan   |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences