

SPOPL Rabbit anti-Human Polyclonal (aa228-357) (FITC) Antibody - LS-C241208 - LSBio	
CatalogID:	LS-C241208
Target:	speckle-type POZ protein-like (SPOPL)
Synonyms:	SPOPL Antibody, Roadkill homolog 2 Antibody, Speckle-type POZ protein-like Antibody, BTBD33 Antibody, HIB homolog 2 Antibody
Host	SPOPL antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	SPOPL antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	SPOPL antibody was raised against kLH-conjugated synthetic peptide (aa228-357) from internal region of human SPOPL.
Specificity:	Human SPOPL
Epitope:	aa228-357
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	